

# VISCO™



## COMING SOON

### Introducing a new type of instrument. Taking you beyond the limits.

- Compact and easily carried with one hand
- Battery operated—measurements can be taken anywhere.
- Fully digital display allows for anyone to easily read results.
- Capable of measurements with only a small amount of sample.
- Simple construction and easy, one-button operation allows for anyone to assemble and operate the unit.
- Any disposable container other than the included beaker can be used.

#### Specifications

|                      |   |
|----------------------|---|
| Cat.No.              | 6800  |
| Measurement Range    | Viscosity A1 50 to 200,000mPa·s, 50 to 200,000cP<br>A2 100 to 600,000mPa·s, 100 to 600,000cP<br>A3 500 to 2,000,000mPa·s, 500 to 2,000,000cP (1mPa·s=1cP) |
| Torque               | 0 to 100% (recommended torque% 10 to 100%)  |
| Temp                 | 10.0 to 40.0°C, 50.0 to 104.0°F   |
| Resolution           | Viscosity more than 10,000 : 1mPa·s, less than 10,000 : 0.1mPa·s<br>Torque 0.1%<br>Temp 0.1°C   |
| Measurement Accuracy | Viscosity ±1%<br>Temp ±0.2°C  |
| Dimensions & Weight  | Main Unit : 12×12×20.06cm 1.2kg<br>Stand : 0.5kg  |

※The above specifications apply only to measurements taken using the standard accessories.

All ATAGO products are designed and manufactured in Japan.

**ATAGO CO., LTD.**

Headquarters: The Front Tower Shiba Koen, 23rd Floor 2-6-3 Shiba-koen, Minato-ku, Tokyo 105-0011, Japan TEL : 81-3-3431-1943 FAX : 81-3-3431-1945  
<http://www.atago.net/> MAIL : [overseas@atago.net](mailto:overseas@atago.net) \* Specifications and appearance are subject to change without notice.

V.01 1507K